




Easy crossing in balanced performance

ACROSS CTO™

OTW

- High pushability
- Seamless transition
- Outstanding crossing
- Flexible balance
- Enhanced trackability




Mid marker

4 mm tip
Seamless Tip-to-shoulder

High RBP 17 Atm

Low entry
profile 0.016"



Ordering Information / Technical Information

ACROSS CTO™ OTW

Reference	Balloon Diameter	Balloon Length	Marker	Tip Length	Entry Profile	Max. Guidewire	Useable Length	Min. Guide Catheter / ID
110100481	1.1 mm	10 mm	1	4 mm	0.016"	0.014"	150 cm	5 F / 0.056"
110150481	1.1 mm	15 mm	1	4 mm	0.016"	0.014"	150 cm	5 F / 0.056"
110200481	1.1 mm	20 mm	1	4 mm	0.016"	0.014"	150 cm	5 F / 0.056"

ACROSS CTO™ flexy

Reference	Balloon Diameter	Balloon Length	Marker	Tip Length	Entry Profile	Max. Guidewire	Useable Length	Min. Guide Catheter / ID
110100482	1.1 mm	10 mm	1	4 mm	0.016"	0.014"	150 cm	5 F / 0.056"
110150482	1.1 mm	15 mm	1	4 mm	0.016"	0.014"	150 cm	5 F / 0.056"
110200482	1.1 mm	20 mm	1	4 mm	0.016"	0.014"	150 cm	5 F / 0.056"

Control your Size

Balloon Ø	Atm	Atm*	Atm	Atm	Atm	Atm	Atm**	Atm
1.1 mm	5	7	9	11	13	15	17	19
1.1 mm	1.04 mm	1.10 mm	1.13 mm	1.16 mm	1.20 mm	1.23 mm	1.28 mm	1.38 mm

Atm* Nominal Pressure (NP) Atm** Rated Burst Pressure (RBP)



The balloon material is a semi-compliant polyamide, providing a 'Work Your Size' growth.

Easy crossing in balanced performance